
Woodrolls Übungen | exercises

Art.nr.: 47020

Art.nr.: 47021


Waden

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calves

rolling back and forth


seitliche

Oberschenkel

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lateral thighs

rolling back and forth


hintere

Oberschenkel

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ham strings

rolling back and forth


innere

Oberschenkel

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inner thighs

rolling back and forth


Gesäß

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bottom

rolling back and forth


vordere

Oberschenkel

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front thighs

rolling back and forth

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
upper back

static holding

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lateral back

rolling back and forth

Schulter

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shoulder

rolling back and forth

Oberarm

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upper arm

rolling back and forth

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
upper back

rolling back and forth

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upper back

rolling left and right

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lower back

rolling back and forth

When exercising for the first couple of times you might find the hardness of the Woodroller a little bit awkward to get accustomed to; this is due to the tackiness and stiffening of the tissue. Things tend to get easier with regular use. It is recommendable to reduce the number of repeats at first and to increase these as time goes on.

The intensity of massage effect can be varied by the pressure exerted on the Woodroller. For example, when exercising the calf or thigh you can lower one leg for support, this catering for ease in performing the exercise. To increase intensity, place one leg over the other. Similarly, you can exert influence on the pressure by shifting the weight of your body or by using the arms and legs for greater support. You will also notice the difference between performing the exercises against a wall and in a reclining posture down on the floor.

Beginners should stick to wall exercises for the time being, gradually stepping things up as more confidence is gained. At first, the muscular system and fascia connective tissue need to get used to the hardness of the material.

Through regular use of the Erzi Woodroller over a longer period of time you will discover the positive impact this has on your physical well-being along with improvement of many functional problems.

On the other hand, if you are suffering from injury or disease of the bones or muscular system, or if for other reasons, there are certain restrictions to your musculoskeletal system, be sure to contact your local doctor or therapist first before attempting these workouts!

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Die **Erzi Woodroll** line enhances fascia training, raising it to a level hitherto unprecedented. Natural beech produces an incredible depth penetration effect on account of its hardness and inherent warmth. In addition, the anti-bacterial properties of natural wood products along with their weight make them ideal for performing soothing stroke-massage applications. Thus all requirements are equally catered for, regardless of whether the rolls are used by beginners or by professional athletes for revitalisation purposes, or as a means of preparing users about to take up some strenuous form of activity. Therapists will simply dote on the high quality standards enhanced by the pleasant feel of wood in its purest state.

These classical fascia rolls cater for enjoyable soothing massage applications as well as for constructive activation and mobility workouts. It is the use of natural beech incorporating inherent warmth along with stability and high quality that makes our fascia rolls so incomparable.

Nothing could be better for sports-event training, daily back-training exercises or just for sheer relaxation brought about by a thorough self-massage. The natural state of the material makes fascia training a very special experience indeed!

Erzi Fascia Woodrolls are the perfect thing for stimulating and relaxing the muscles before engaging in any strenuous activity or purely and simply for revitalisation. On account of their large-dimensioned bearing surface and the gentle pressure exerted, the rolls are ideal for users preparing for other Erzi Woodroll shapes producing a greater depth penetration effect.

When exercising, always keep the abdomen, back and backside in a tensed state. When performing exercises in a reclining posture down on the floor, avoid rolling too far away from the body's centre of gravity so as not to lose your balance. We recommend 8 to 12 repeats per exercise. However, this is merely a suggestion. You can perform as many exercises as you wish.

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**Erzi®** Qualitätsprodukte aus Holz GmbH